



MEGHdin


Povidone Iodine

Non-adherent wound dressing

10% U.S.P. (available Iodine 1%)

- ✓ Easy to use
- ✓ Antimicrobial action
- ✓ Cost effective





Non-adherent Dressing

- Enhances patient comfort
- Minimises potential damage to granulation tissue

Antimicrobial Action

- Broad spectrum antimicrobial action of Povidone-Iodine (PVP-I)
- Effective against a wide range of pathogens¹


Colour indicator for dressing change

INTENDED USE


- ✓ Minor burns
- ✓ Skin tears
- ✓ Diabetic foot ulcers
- ✓ Minor traumatic skin loss injuries
- ✓ Ulcers deriving from different aetiologies
- ✓ Pressure ulcers

INFECTION PREVENTION

MEGHdin offers comprehensive antibacterial efficacy against a range of different pathogens, including:




Bacteria, including MRSA



Fungi



Yeasts



Protozoa



Product Size	Units Per Box	MPC
5cm x 5cm	25	P05
9.5cm x 9.5cm	10	P95

1. Data on file